Test destiné à l'autodiagnostic.

INTRODUCTION

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Test to be used for self-diagnosis.

INTRODUCTION

COVID-19 is an acute respiratory disease, caused by infection via the SARS-CoV-2 virus. The SARS-CoV-2 virus belongs to the same family as the viruses which cause Severe Acute Respiratory Syndromes (SARS).

The most common symptoms of COVID-19 are fever, fatigue and a dry cough. Certain patients may experience generalised pain, nasal congestion, rhinorrhoea (runny nose), a sore throat, diarrhoea, and loss of taste and smell. These symptoms are typically mild, and develop gradually. Some people may be infected by the virus, but present no symptoms. The majority of people (approx. 80%) recover from the virus without necessitating hospitalisation.

The disease can be passed on from person to person by tiny droplets, which are transmitted from the nose or mouth when someone infected with SARS-CoV-2 coughs or simply exhales. Transmission by contact of contaminated hands touching the face is also frequent. Effectively, hands may be contaminated with the virus following contact with surfaces carrying the virus. The estimated period of incubation for the SARS-CoV-2 virus varies from one (1) to fourteen (14) days.

INTENDED USE

The Biosynex COVID-19 Ag Antigenic Self-testing Kit is a rapid, in vitro immunochromatography assay (also known as a lateral flow test) which allows for the qualitative detection of antigens of the nucleocapsid protein (N) of the SARS-CoV-2 virus, in samples of nasal swabs. It has been designed to assist in rapid diagnosis of infections of the SARS-CoV-2 virus. This test is a qualitative, rapid immunochromatography assay, which uses highly sensitive monoclonal antibodies to detect the nucleocapsid protein (N) of the SARS-CoV-2 virus.

This self-testing kit allows you to identify whether you are, in fact, infected by the virus which causes COVID-19. This test is only to be used by adults, or under the supervision of an adult if the test is to be carried out on a child.

A negative result from the Biosynex COVID-19 Antigenic Self-Testing Kit for an asymptomatic individual does not wholly rule out a COVID-19 infection. If the individual presents symptoms suggestive of a COVID-19 infection, a negative result must be verified with a laboratory test.

KIT COMPOSITION

Ref 859257 contains one bag for one self-test.
Ref 859262 contains five bags, one self-test can be run with each bag.

SAMPLING ADVICE

- Wash hands before each test, in order to reduce the risk of contamination.
- For accurate results, do not use any samples which appear too viscous, or which visibly contain blood. You can blow your nose to remove any excess mucus before completing the test.

TEST LIMITS

- This test does not determine the aetiology of the respiratory infection caused by micro-organisms other than the SARS-CoV-2 virus. The Biosynex COVID-19 Ag Antigenic Self-testing Kit can detect both the viable and the non-viable SARS-CoV-2 virus.
- The test must be used to detect the antigen of the SARS-CoV-2 virus, using a nasal swab. Neither the quantitative value nor the rate of increase for the concentration of the SARS-CoV-2 virus can be determined using this qualitative test.
- The accuracy of the test depends on the quality of the swab sample. False negative results may be given following poor sampling.
- Any failure to respect the test procedure may negatively impact the performance of the test and/or invalidate the test result.
- If the result of the test is negative, yet clinical symptoms persist, it is advised that you carry out additional tests using other clinical methods. A negative result at no time rules out the presence of antigens of the SARS-CoV-2 virus in the sample, as they may be present but at a level inferior to the minimum detection level of the test, or if the sample has been collected incorrectly.
- A negative result does not rule out infection by the SARS-CoV-2 virus, particularly in people who have come into contact with the virus. Follow-up tests with molecular diagnostics should be scheduled to rule out infection in these people.
- This test is not a substitute for a medical consultation, or for the result of a biological analysis carried out in a medical analysis laboratory.
- Positive test results do not exclude the possibility of co-infections of other pathogens.
- Positive test results do not allow for differentiating between the SARSCoV virus and the SARS-CoV-2 virus.
- The antigen detected by the test is the N protein. The different variants of the virus reported, at the time of publishing, in certain countries (United Kingdom, South Africa, Brazil, etc.) concern mutations of the Spike protein and therefore have no impact on the functionality of the test.
- Erroneous errors may be obtained:
 - If the test is not used in line with the instructions provided,
 - If the aluminium sachet is damaged, or if the test has not been carried out immediately after opening the aluminium sachet,
 - If the storage conditions are not respected, or if the test is carried out after the expiry date provided on the aluminium sachet.
- A positive result must be confirmed by a laboratory analysis. Consult your doctor, and do not make any medical decisions without first seeking advice from your GP.

HANDLING INSTRUCTIONS

- Keep the self-testing kit and its components out of reach of children; the diluent may pose a significant risk if ingested.
- The used test cassette, as well as all components of the test, can be disposed of along with household waste, in a securely closed bag.
- All components included in this kit are intended solely for use for this test. Do not re-use this test and/or any of its components.
- This test must be carried out at room temperature (15 °C – 30 °C) and in an environment free from excessive humidity.
- This is a single-use test.
- This test must be used immediately after collecting the sample.
- The extraction diluent contains a solution containing a preservative (0.09% sodium azide). If the solution comes into contact with skin or eyes, rinse thoroughly with water.
- Solutions containing sodium azide may react explosively with lead or copper pipes. Use large volumes of water to thoroughly rinse waste disposed of in a sink.
- Do not use the cassette if it has been removed from its sealed packaging for over an hour.
- Do not use the cassette if its packaging is damaged.

STORAGE

- This test must not be frozen, and must be stored in a dry place between 2 °C and 30 °C.
- The test cassette must stay in the sealed sachet until it is about to be used.
- Do not use the self-testing kit after the expiry date provided on the aluminium sachet, and on the outside of the cardboard box.

SAMPLE COLLECTION

Scan the QR code below to head over to YouTube to watch the demonstration video.

Please read the notice in its entirety. Have a watch or timer to hand.

- Wash your hands before handling the kit to avoid any potential contamination.
- Open the box and remove **all elements**. Place them on a clean, flat surface. Identify each element of the kit.
- Tear off the top of the pouch, following the indent. Remove the cassette. **Use the test within the hour.**

- Remove the seal of the tube containing the diluent, and place it in the cut-out hole in the box.
- Remove the swab from its packaging.

Samples to be taken from each nostril using the same swab:

- Insert the swab vertically into a nostril, **until it starts to resist** (approx. 3 cm).
- Roll the swab **5 times** along the mucus inside the nostril, to ensure that the mucus and cells are collected. **The swab must touch the inner walls of the nostrils.**
- Using the same swab, **repeat this process in the other nostril**. Remove the swab from the nasal cavity.

- Insert the swab into the tube and rotate it at least **6 times** whilst pressing the end of the swab against the bottom of the extraction tube.
- Leave the swab in the extraction tube **for 1 minute**.
- Extract as much liquid as possible from the swab by "squeezing it out" using the edge of tube.** To do this, press the swab against the edge of the tube several times. Then, remove the swab from the tube and dispose of it.


COMPLETION

- Attach the dropper teat onto the tube.
- Using the dropper, **slowly add 4 drops** of solution into well S.
- Start the timer, or note down the time of deposit. **Read the results after 15 minutes.** Do not read after 20 minutes.
- Place the red stopper over the dropper teat before discarding the tube into the bin, in order to avoid any potential contamination.

Refer to the "Interpretation" section for details on test reading.

Accessory	Manufacturer	Proxy	CE Marking
Swab	Jiangsu Changfeng Medical Industry Co., Ltd Touqiao Town, Guangling District, Yangzhou, Jiangsu 225109 CHINA	Llins Service & Consulting GmbH Obere Seegasse 34/2, 69124, Heidelberg, GERMANY	CE0197

INTERPRETATION

The test result is **NEGATIVE** if a single coloured band only appears in the control zone (C), and no coloured band appears in the test zone (T), **this applies regardless of the intensity of the band**.
A negative result indicates the absence of antigens of the SARS-CoV-2 virus.
A negative result does not rule out recent infection of the SARS-CoV-2 virus. If you think that you have come into contact with the virus (with an infected person) over the days preceding the completion of the test, we would advise you to complete another test with a laboratory. Continue to protect yourself and maintain social distancing.

The result is **POSITIVE** if 2 coloured bands appear on the membrane. One coloured band will appear in the control zone (C), and another will appear in the test zone (T), **this applies regardless of the intensity of the bands**. A positive result indicates that you are currently carrying the SARS-CoV-2 virus.
Consult your doctor, continue to protect yourself, and maintain social distancing.

If the control line (C) does not appear after having carried out the test, the result is considered to be **INVALID**.
Review the procedure and repeat the test with a new cassette.
An insufficient sample volume or an excessively viscous sample, incorrect procedural techniques (incorrect swab collecting, temperature and humidity conditions in which the test was carried out), or tests which have been opened for more than an hour, or which have expired, are the most common reasons for the lack of control line appearance.

To discard of your testing kit and its components, please refer to the guidelines in force in your local area.

PERFORMANCE


Sensitivity and specificity
The Biosynex COVID-19 Ag Antigenic Self-testing Kit has been assessed using patient samples. A marketed molecular test was used as the reference method. The study included 249 samples (109 confirmed positive, and 140 negative samples).

	%
Sensitivity	97.2%
Specificity	100 %

A feasibility study demonstrated that:
- 97,0% of different types of results have been correctly interpreted.
- 87,0% of laypersons have completed the test without need for assistance.

INTERFERENCES

No positive or negative interference has been exhibited by the study in a matrix of nasopharyngeal secretions with the following substances: human blood (with EDTA anticoagulant), Mucin, antiviral medication (Osetamivir phosphate, ribavirine), Antibiotics (Levofloxacin, Azithromycin, Meropenem, Tobramycin), nasal sprays or drops (Phenylephrine, Oxymetazoline, Alkalol nasal wash, 0.9% NaCl), nasal corticosteroids (Beclomethasone, Hexadecadrol, Flunisolide, Triamcinolone, Budesonide, Mometasone, Fluticasone, propionate of fluticasone).

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