

## Rapid

Results in minutes

## Convenient

Easy Handling

## Reliable

Constant Quality

## Efficient

Competitive

# IMMUNOQUICK®

## Adenovirus, Rotavirus and Combo

The Rota- and Adenoviruses are both RNA viruses. These viruses account for most cases of acute and severe diarrhea in young children.

The Rotavirus, with an incubation period of between 1-3 days, is highly contagious and transmitted through the fecal-oral route.


Besides diarrhea the patient infected with the Adenovirus may have fever and may be vomiting. Other symptoms of an infection may include, but are not limited to: conjunctivitis, cystitis or respiratory diseases.

Rota- and Adenovirus infections are mostly self-limited illnesses for persons

with healthy immune systems. The treatment, in case of gastroenteritis in conjunction with diarrhea, consists of increased fluid intake to avoid dehydration.

The IMMUNOQUICK® Adeno-, Rotavirus and Combo are lateral flow immunoassays for the fast and easy detection of these viruses in stool samples. These tests are used to differentiate Rota- or Adenovirus infections from other pathogens causing similar symptoms, which is an aid to the physician to decide on the right treatment.

## Interpretation



### POSITIVE RESULT

When in the test (T) zone a red and the control (C) zone a blue band appears.



### NEGATIVE RESULT

Only one blue band appears in the control (C) zone .



### INVALID RESULTS

When there is no band visible in the control (C) zone the test is invalid and should be repeated with a new device.



## Additional Information

**Guaranteed shelf life:** 12 months

**Storage conditions:** 2-30°C

### IMMUNOQUICK Adenovirus

**REF 0540\_K25:** each kit contains 25 tests

### IMMUNOQUICK Rotavirus

**REF 0541\_K25:** each kit contains 25 tests

### IMMUNOQUICK Adeno/Rota Combo

**REF 0542\_K25:** each kit contains 25 tests

## Performance

	Sensitivity	Specificity
IMMUNOQUICK® Adenovirus	>99,9%	99,4%
IMMUNOQUICK® Rotavirus	99,1%	>99,9%

**Biosynex**, together with **DIMA GmbH**, is a member of the Biosynex group. The Biosynex group is one of the most significant manufacturers of rapid membrane tests on the in-vitro diagnostic market.

We offer high quality, easy to use Drugs of Abuse, Tropical Disease, Virology, Bacteriology, Autoimmune Disease, Oncology, Cardiology, Fertility, Obstetrics and Emergency Medicine kits.

Our Research and Development laboratories in Germany and France and our state of the art manufacturing facilities are dedicated to continuously meet the needs of our customers.

Our aim is to provide diagnostic tests to help physicians make faster and better medical and therefore lifesaving decisions.

BIO SYNEX